

## Checklist of the practical skill "Applying hemostatic twists"

FULL NAME. student \_\_\_\_\_ Group, faculty \_\_\_\_\_ Date \_\_\_\_\_

№	Act	1*	0.5**	0***
1	Examine the area of damage and check for arterial bleeding	5.0	3.0	0
2	Choose the right place for the application of the twist proximal to the wound and as close as possible	5.0	3.0	0
3	Identify contraindications for applying a twist	5.0	3.0	0
4	Raise limb	10.0	7.0	0
5	Place a soft cloth (bandage, clothing, etc.) on the place where the curl is applied,	10.0	4.0	0
6	Correct swirling between two nodes. The twist is tightened so that the bleeding stops, but not overtighten	25.0	15.0	0
7	The stick is securely fixed	15.0	10.0	0
8	After applying the twist, the bleeding stopped, there is no pulsation in the periphery	15.0	-	0
9	A note is attached to the twist, indicating the time of its imposition	10.0	5.0	0
	Total	100.0	50.0	0

\* completed completely; \*\* incompletely executed (with errors); \*\*\* did not comply

